

Safety Guidelines

Chewable Products and Chewies

⚠ CAUTION All Fun and Function chewable products are for therapeutic use only and they should not be considered toys. Please read these Safety Guidelines and the product insert before using your chewy. Consult with your pediatrician or therapist for oral motor guidance.

USING YOUR CHEWY

- Adult supervision is always required when using chewable products and chewies.
- Caregiver should have knowledge of each child's unique needs and developmental capabilities. Some children require additional supervision for their own safety.
- Caregiver should teach child how to correctly use the chewy:
 - Only mouth the chewable portion.
 - Do not chew on necklace cords or clasp.
- Necklace cords have a safety clasp to prevent strangulation and to breakaway when pulled on. Never use if safety clasp is broken.
- Remove chewy before active play, naptime, and bedtime.

INSPECTING YOUR CHEWY

- No product is indestructible. Before each use, inspect chewy and cord for signs of wear and tear.
 - If the user has aggressive biting or chewing habits, inspect more frequently.
- Immediately discard and replace
 - Any chewy that shows signs of discoloration, cracking, or any other signs of deterioration or damage.
 - Any cord that shows signs of fraying, thinning, breakage or damage.

CARING FOR YOUR CHEWY

- After use, hand wash or surface wash the chewable product according to the care instructions provided. Air dry.
- Some chewable products are dishwasher safe, top rack. Check specific product information or email us at wecare@funandfunction.com.

ABOUT OUR CHEWIES

- Our chewy products have been tested to withstand varying levels of pressure, for at least 2000 repetitive cycles.
- Chewable product durability may depend on the oral sensory needs of the child, as well as oral habits. If the user has aggressive biting or chewing habits, choose a higher durability chewy and inspect more frequently.

At Fun and Function, safety is our priority. If you have questions about how to safely use this product or any Fun and Function chewable product, please contact us at



fun & function

www.funandfunction.com/safetyguidelines

Fun and Function LLC, PO Box 11, Merion Station, PA 19066

Safety Guidelines

Sensory Room Products

(including SensaSoft™, Bubble Tubes, Tactile Wall Panels and Fiber Optics)

⚠ CAUTION All Fun and Function sensory products are for therapeutic use only and they should not be considered toys. Equipment is highly specialized and is not a toy. Please read these Safety Guidelines and the product insert before using your sensory room products.

INSTALLING YOUR SENSORY ROOM PRODUCTS

Read the product insert thoroughly for proper installation instructions.

- Sensory Room products are for indoor use only, unless otherwise specified.
- Permanent installations should be done by a professional installer.
 - A professional is a certified construction professional, in compliance with applicable regulations in the province, state, or country where the product is installed. We recommend that you find a professional with appropriate licensing and proof of liability insurance coverage. Installer must follow standard safety practices for installation as outlined in the product insert.
- To prevent tip over and fall hazards, appropriate installation hardware must be used including required brackets and anchors.
- Ensure any electrical cords and outlets are secure and inaccessible to users.

USING YOUR SENSORY ROOM PRODUCTS

Adult supervision is always required when using sensory room products.

- Consult with your pediatrician or therapist to confirm that the user's developmental abilities are correctly matched to the equipment.
- Tip over and fall hazard: ensure all equipment is properly anchored.

INSPECTING AND MAINTAINING YOUR SENSORY ROOM PRODUCTS

- Inspect equipment prior to each use.
- Carry out a routine safety check at regular intervals to inspect all components including padded structures, LED components, fiber optic strands and liquid filled chambers.
- Discontinue use immediately if product exhibits wear and tear or damage.
- Keep product inserts in an accessible place so all users are familiar with proper use and inspection requirements.

ABOUT OUR SENSORY ROOM PRODUCTS

- Sensory Room products meet all applicable safety requirements for flammability*, banned phthalates, lead, and physical and mechanical hazards. Electronic components are UL certified for safety. (*Check your jurisdiction's building code requirements.)

At Fun and Function, safety is our priority. If you have questions about how to safely use this product or in choosing appropriate sensory room products, please contact us at wecare@funandfunction.com



fun & function

www.funandfunction.com/safetyguidelines

Fun and Function LLC, PO Box 11, Merion Station, PA 19066

Safety Guidelines

Swing and Suspension Products

⚠ CAUTION All Fun and Function products are for therapeutic use only and they should not be considered toys. Please read these Safety Guidelines and the product insert before using your swing or suspension product. Consult with your physician or therapist for motor skills guidance or if user is ill, injured, has respiratory difficulty, or has unique physical or developmental needs.

USING YOUR SWING OR SUSPENSION PRODUCT

Active eyes-on adult supervision is required when using all swings and suspension products.

- **⚠ WARNING** Fall and strangulation hazards may occur if products are used incorrectly.
- Most hammock swings are for gentle rocking and not active swinging.
- Remove or secure swing away from children after each use. Children may become entangled and injured in swing ropes, soft fabric, and hammock swings.
- Ensure that child is an appropriate weight range for the product's intended use.
- Ensure that protective floor mats of appropriate thickness are installed under and 6 feet from the perimeter of the equipment.

INSTALLING YOUR SWING OR SUSPENSION PRODUCT

Read the product insert thoroughly for proper installation instructions

- Swings and suspension products are for supervised indoor use only unless specified. Choose a location where there is a minimum clearance of 6 feet around the perimeter of the swing, unless other distance is specified.
- For doorway bar installations that do not permit a 6-foot clearance, only 2-point suspension products can be used.
- Permanent installations should be done by a professional installer.
 - A professional is a certified construction professional, in compliance with applicable regulations in the province, state, or country where the product is installed. We recommend that you find a professional with appropriate licensing and proof of liability insurance coverage. Installer must follow standard safety practices for installation as outlined in the product insert.
- Ensure appropriate hardware is used including Fun and Function rotational devices, swivels or two-to-one converters. For single point suspension, use with swivel. Hardware requirements are specified in each product listing and instructions.
- Extend protective floor mats of appropriate thickness under and 6 feet from the perimeter of the equipment.

INSPECTING AND MAINTAINING YOUR SWING OR SUSPENSION PRODUCT

- Inspect equipment prior to each use.
- Carry out a routine safety check at regular intervals to inspect all components including fabric, rope, carabiners, and attachment points on both the product and building structure.
- Discontinue use immediately if product exhibits wear and tear or damage.
 - Regular wear and tear is expected of all components, including but not limited to eyebolts, carabiners, O-rings, hooks, swivels, chains and ropes.
 - Replacement parts can be purchased directly from Fun and Function.
- Keep product inserts in an accessible place so all users are familiar with proper use and inspection requirements.

ABOUT OUR SWINGS AND SUSPENSION PRODUCTS

- Swing and suspension products have been tested for physical and mechanical hazards, load bearing and weight capacity. Weight capacity differs depending on product and installation method.

At Fun and Function, safety is our priority. If you have questions about how to safely use this product or in choosing an appropriate swing, suspension products or associated hardware, please contact us at wecare@funandfunction.com.



fun & function

www.funandfunction.com/safetyguidelines

Fun and Function LLC, PO Box 11, Merion Station, PA 19066

Safety Guidelines

Gel-Filled Products

⚠ CAUTION All Fun and Function gel-filled products are for therapeutic use only and they should not be considered toys. Please read these Safety Guidelines and the product insert before using your gel-filled products. Consult with your pediatrician or therapist for developmental and motor skills guidance.

USING YOUR GEL-FILLED PRODUCTS

Adult supervision is required when using all gel-filled products including hand-held, floor and seat products.

- Gel-filled products can become a hazard if product is pierced or damaged.
- Gel-filled products should never be chewed on, ingested, or swallowed.
- Follow the age guidelines stated on the product insert.

INSPECTING GEL-FILLED PRODUCTS

- Inspect gel-filled products prior to each use.
- Discard gel-filled products that show tears, gel leaks, wear at edges and punctures.

CARING FOR GEL-FILLED PRODUCTS

- Never heat or cool any gel-filled product unless specified in instructions. Temperature extremes can compromise product safety.
- Gel-filled products are not designed to be thrown or dropped with excessive force.
- Keep sharp objects away from product. Do not puncture or cut into gel-filled products
- Surface wash and air dry. Products can be washed with mild disinfectant if required.
- Keep product inserts in an accessible place so all users are familiar with proper use and inspection requirements.

ABOUT OUR GEL-FILLED PRODUCTS

- This gel-filled product, including the outer and inner materials, meets all applicable safety requirements for lead, phthalates, physical and mechanical properties, and flammability.
- Gel is non-toxic. User should drink water if swallowed. If gel gets in eyes, flush eyes with water. As when any non-food item is swallowed or gets in eyes, call a doctor.

At Fun and Function, safety is our priority. If you have questions about how to safely use this product or in choosing appropriate gel-filled products, please contact us at wecare@funandfunction.com.



fun & function

www.funandfunction.com/safetyguidelines

Fun and Function LLC, PO Box 11, Merion Station, PA 19066

Safety Guidelines

Putty Products

⚠ CAUTION All Fun and Function putty products are for therapeutic use only and they should not be considered toys. Please read these Safety Guidelines and the product insert before using your putty products. Consult with your pediatrician or therapist for developmental and motor skills guidance.

USING YOUR PUTTY PRODUCTS

Provide active, eyes-on adult supervision so child avoids putting putty in ear, hair, on furniture, appliances, clothing or carpeting.

- Putty products containing small parts can become a choking hazard. Putty products are NOT recommended for children under 3 years of age. Follow the age guidelines stated on the product insert.
- Putty products are not intended for mouthing or chewing and should NEVER be ingested or swallowed.
- Store putty products away from children after each use.

INSPECTING PUTTY PRODUCTS

- Inspect putty products prior to each use.
- Look for wear and tear, cleanliness or material breakdown. Replace as needed.

CARING FOR PUTTY PRODUCTS

- Never heat any putty product. Temperature extremes can compromise product safety.
- Keep product inserts in an accessible place so all users are familiar with proper use and inspection requirements.

ABOUT OUR PUTTY PRODUCTS

- Putty products come in different resistance levels: extra soft/soft/medium/firm. Resistance level reflects user's specific needs such as muscle tone or sensory needs. Consult your therapist for guidance.
- This putty product meets all applicable safety requirements for lead and physical and mechanical hazards.

At Fun and Function, safety is our priority. If you have questions about how to safely use this product or in choosing appropriate putty products, please contact us at wecare@funandfunction.com.



fun & function

www.funandfunction.com/safetyguidelines

Fun and Function LLC, PO Box 11, Merion Station, PA 19066

Safety Guidelines

Seated Scooters

⚠ CAUTION All Fun and Function products are for therapeutic use only and they should not be considered toys. Please read these Safety Guidelines and the product insert before using your scooter. Consult with your physician or therapist for motor skills guidance or if user is ill, injured, has respiratory difficulty, or has unique physical or developmental needs.

USING YOUR SEATED SCOOTER

Active eyes-on adult supervision is required when using all seated scooters.

- Allow only one rider at a time.
- Do not exceed the weight capacity indicated on the product insert.
- Fall hazard: Do not allow child to stand on the seated scooter. Never allow child to use seated scooters on or near stairs or other hazards.
- Store seated scooters away from children after each use.

INSPECTING AND MAINTAINING YOUR SCOOTER

- Some scooters require assembly. Follow directions as outlined in the product insert.
- Inspect equipment prior to each use for cracks, breaks or irregular wear.
- Carry out a routine safety check at regular intervals to inspect all components including wheels/casters and body of scooter.
- Discontinue use immediately if product exhibits wear and tear or damage.
 - Some regular wear and tear is expected. Replacement wheels can be purchased directly from Fun and Function.
- Scooters are for indoor use only.
- Surface wash, air dry.
- Keep product inserts in an accessible place so all users are familiar with proper use and inspection requirements.

ABOUT OUR SCOOTERS

- Scooters meet all applicable safety requirements for physical and mechanical hazards, lead, and phthalates.
- Weight capacity and recommended age range differs depending on product. See product insert for details.

At Fun and Function, safety is our priority. If you have questions about how to safely use this product or in choosing an appropriate scooter, please contact us at wecare@funandfunction.com.



fun & function

www.funandfunction.com/safetyguidelines

Fun and Function LLC, PO Box 11, Merion Station, PA 19066

Safety Guidelines

Weighted Clothing & Products (excluding blankets)

(including vests, lap pads, belts, bean bags, hats and other therapeutic items)

⚠ CAUTION All Fun and Function weighted products are for therapeutic use only and they should not be considered toys. Please read these Safety Guidelines and the product insert before using your weighted products. Consult with your pediatrician or therapist for developmental and motor skills guidance.

USING YOUR WEIGHTED PRODUCTS

- Adult guidance is recommended when using all weighted clothing and products. Stop use if user reports fatigue, pain, or any discomfort.
- Follow the age and weight guidelines
 - User should be at least 3 years of age, unless age guidelines state otherwise.
 - Weight should be adjusted for each user's needs.
 - * 10% of total body weight is a recommended standard.
 - * Total weight should never exceed 10% of total body weight.
- Weighted products should be incorporated into use gradually.
 - Begin by wearing the product with lower weights and for short periods of time (e.g. 10 - 20 minutes) and then removing the weight.
 - The amount of time the product is worn and the amount of weight added to the product can be increased gradually, as directed by a professional, such as an occupational therapist.
- Only one weighted product should be used at a time.
- Do not allow user to sleep in weighted clothing. Only weighted blankets are suitable for sleep.
- User should be able to remove weighted product without help.
- Consult physician before using if user is ill, injured, has respiratory difficulty, or has unique physical or developmental needs.
- Never use to restrict movement, restrain, or cover user's face.

INSPECTING YOUR WEIGHTED PRODUCTS

- Inspect weighted products prior to each use. Inspect both the outer product and inner weighted portion if removable.
- **⚠ WARNING:** Choking hazard. Look for loose stitching, ripped seams, fraying, torn or worn areas that may release weighted filling.
- If product is damaged, immediately stop use until it is properly repaired or replaced.

CARING FOR YOUR WEIGHTED PRODUCTS

- Follow cleaning instructions as outlined on the product insert. Some products can be machine or handwashed, others are surface clean only. Air dry.
- Keep product inserts in an accessible place so all users are familiar with proper use and inspection requirements.

ABOUT OUR WEIGHTED PRODUCTS

- This weighted product has been tested for lead, phthalates, physical and mechanical properties, cleanliness, and flammability.
- The product is designed and manufactured for long-term durability for our most active users.

At Fun and Function, safety is our priority. If you have questions about safe play, safe practices, or choosing appropriate weighted products, please contact us at wecare@funandfunction.com.



fun & function

www.funandfunction.com/safetyguidelines

Fun and Function LLC, PO Box 11, Merion Station, PA 19066

Safety Guidelines

Weighted Blankets

⚠️ **CAUTION** All Fun and Function weighted blankets are for therapeutic use only. Consult with your pediatrician or therapist for developmental and motor skills guidance.

USING YOUR WEIGHTED BLANKET

Adult guidance is recommended when using all weighted blankets. Stop use if user reports fatigue, pain or any discomfort.

- Follow the age and weight guidelines
 - User should be at least 3 years of age.
 - Weight should be adjusted for each user's needs.
 - * 10% of total body weight is a recommended standard.
 - * Total weight should never exceed 10% of total body weight.
- User should be able to remove blanket without help.
- Consult physician before using if user is ill, injured, has respiratory difficulty, or has unique physical or developmental needs.
- Never use more than one weighted product at a time.
- Never use to restrict movement, restrain, or cover user's face.
- Never use in a crib.

INSPECTING YOUR WEIGHTED BLANKET

- Inspect weighted blanket prior to each use.
- ⚠️ **WARNING: Choking Hazard:** Look for loose stitching, ripped seams, fraying, torn or worn areas that may release weighted filling.
- If blanket is damaged, immediately stop use until it is properly repaired or replaced.

CARING FOR YOUR WEIGHTED BLANKET

- Machine or hand wash in cold water depending on product recommendation. Air dry.
- Keep product inserts in an accessible place so all users are familiar with proper use and inspection requirements.

ABOUT OUR WEIGHTED BLANKET

- This weighted blanket has been tested for lead, phthalates, physical and mechanical properties, cleanliness, and flammability.
- The product is designed and manufactured for long-term durability for our most active users.

At Fun and Function, safety is our priority. If you have questions about safe play, safe practices, or choosing appropriate weighted products, please contact us at wecare@funandfunction.com.



fun & function

www.funandfunction.com/safetyguidelines

Fun and Function LLC, PO Box 11, Merion Station, PA 19066

Safety Guidelines

Stretch Fabric Products

(including spandex therapy tools, transformer sensory sack, hug sleeping bag and other therapeutic items; excluding stretch clothing)

⚠️ CAUTION All Fun and Function stretch fabric products are for therapeutic use only and they should not be considered toys. Please read these Safety Guidelines and the product insert before using your stretch fabric products. Consult with your pediatrician or therapist for developmental and motor skills guidance.

USING YOUR STRETCH FABRIC PRODUCTS

Adult supervision is required when using all stretch fabric products.

- ⚠️ WARNING Fall Hazard:
 - Never suspend a person in a stretch product/spandex sensory tool or wrap around feet.
 - Use extra caution on uncovered floors. Smooth material is slippery.
 - Use appropriate spotting techniques when a person is in the spandex sensory tool.
- ⚠️ Suffocation Hazard: Do not cover head for prolonged periods of use.
- Only compression snuggle sheets and sleeping bags are suitable for sleep. Do not use other stretch fabric products for sleep.
- Follow the age and size guidelines stated on the product insert.

INSPECTING YOUR STRETCH FABRIC PRODUCTS

- Inspect stretch fabric products prior to each use.
- Look for loose stitching, ripped seams, fraying, torn or worn areas.
- If product is damaged, immediately stop use until it is properly repaired or replaced.

CARING FOR YOUR STRETCH FABRIC PRODUCTS

- Handwash and air dry unless otherwise specified.
- Keep product inserts in an accessible place so all users are familiar with proper use and inspection requirements.

ABOUT OUR STRETCH FABRIC PRODUCTS

- Stretch fabric products have been designed to support children with various sensory needs.
- This stretch fabric product meets all applicable safety requirements for physical and mechanical properties, flammability and seam strength.
- The product is designed and manufactured for long-term durability for our most active users.

At Fun and Function, safety is our priority. If you have questions about how to safely use this product or in choosing an appropriate stretch fabric product, please contact us at wecare@funandfunction.com.



fun & function

www.funandfunction.com/safetyguidelines

Fun and Function LLC, PO Box 11, Merion Station, PA 19066

Safety Guidelines

Bean Bag and Foam Filled Products

(including Bean Bag Chairs, Crash Pads, Mushy Smushy and mini Mushy Smushy products)

⚠ CAUTION All Fun and Function bean bag products are for therapeutic use only and they should not be considered toys. Please read these Safety Guidelines and the product insert before using your bean bag products. Consult with your pediatrician or therapist for developmental and motor skills guidance.

USING YOUR BEAN BAG CHAIR/MUSHY SMUSHY

Adult supervision is required when using all bean bag products.

- ⚠ WARNING Suffocation and Choking Hazard:
 - Children can suffocate if face is pressed into beanbag.
 - Children can suffocate or choke on bean bag filling if accessed.
 - Children under 2 should never be placed on or around bean bag products.
 - Bean bag products should never be used for sleeping.
- Children should never be allowed to play underneath bean bag products or climb inside bean bag product covers.
- Follow the age and size guidelines stated on the product insert - most are for ages 3 and older.

INSPECTING BEAN BAG PRODUCTS

- Inspect bean bags prior to each use.
- Look for exposed filling, loose stitching, ripped seams, fraying, torn or worn areas that may release filling.. If product has zipper cover, ensure zipper is working functionally.
- If product is damaged, do not use until it is repaired or replaced.

CARING FOR YOUR BEAN BAG PRODUCTS

- Spot wash and air dry.
- If product has removable cover, cover can be washed and air dried. Cover can only be removed using a tool, such as a paperclip. No zipper pull is included to prevent children from crawling inside of beanbag. After washing, ensure cover and zipper are firmly reaffixed so child cannot access inside of beanbag.
- Keep product inserts in an accessible place so all users are familiar with proper use and inspection requirements.
- Keep sharp objects away from bean bag to avoid puncturing.

ABOUT OUR BEAN BAG PRODUCTS

- Bean Bag products have been designed to support children with various sensory needs.
- This bean bag product has been tested for physical and mechanical properties, flammability and seam strength.
- The product is designed and manufactured for long-term durability for our most active users.

At Fun and Function, our family's safety is our priority. If you have questions about safe play, safe practices, or choosing appropriate bean bag products, please contact us at wecare@funandfunction.com.



fun & function

www.funandfunction.com/safetyguidelines

Fun and Function LLC, PO Box 11, Merion Station, PA 19066

Safety Guidelines

Climbing Wall Products

⚠ CAUTION All Fun and Function products are for therapeutic use only and they should not be considered toys. Please read these Safety Guidelines and the product insert before using your climbing wall product. Consult with your physician or therapist for motor skills guidance or if user is ill, injured, has respiratory difficulty, or has unique physical or developmental needs.

USING CLIMBING WALL PRODUCT

- **⚠ WARNING** Fall Hazard: Serious personal injury, disfigurement, or death may occur due to the inherent risks associated with climbing.
- Active eyes-on adult supervision is required when using all climbing wall products.
- Choose products and hardware with an appropriate weight range for the product's intended use.
- Close climbing wall areas from child access when supervised use is concluded.

INSTALLING YOUR CLIMBING WALL PRODUCT

- With any of our products which require installation and which bear full or partial weight of the user (e.g., swings, climbing walls, etc.), we strongly recommend that these products be installed by a qualified individual and/or by a professional installer. If the products detach from a ceiling and wall, death or serious injury can occur. For example, it is possible that an installer may not properly use mounting screws affixed to the ceiling beam or wall structure (wood studs, metal studs, or masonry) or that "hollow wall" or other anchors may be incorrectly used. These installations could fail.
- To prevent this hazard, a qualified person must install the products, as indicated in the instructions.
 - A qualified individual has a good understanding of the ceiling or wall structure where the product is to be anchored, including the positioning of the wall studs and structural components in the case of a wood or metal structure.
- If the structure of the ceiling or wall is not a wood stud, metal or masonry or if the installer has any doubt about how or where to anchor the products, then the service of a professional is absolutely required for anchoring the products.
 - A professional is a certified construction professional, in compliance with applicable regulations in the province, state, or country where the product is installed. We recommend that you find a professional with appropriate licensing and proof of liability insurance coverage.
- If you have any doubt about installation requirements, do not attempt to install the product. Please contact us with requests for returns.

INSPECTING AND MAINTAINING YOUR CLIMBING WALL PRODUCT

- Inspect equipment – including the safe placement of fall mats – prior to each use.
- Carry out a routine safety check at regular intervals to inspect all components including climbing holds and attachment points on both the product and building structure.
- Discontinue use immediately if product exhibits wear and tear or damage.
 - Some amount of ordinary wear and tear is expected of all components, including but not limited to climbing holds and screws. Replacement parts can be purchased directly from Fun and Function.
- Keep product inserts in an accessible place so all users are familiar with proper use and inspection requirements.

ABOUT OUR CLIMBING WALL PRODUCTS

Climbing wall products have been tested for physical durability, mechanical hazards, and weight capacity. Weight capacity differs depending on product and installation method.

At Fun and Function, safety is our priority. If you have questions about how to safely use this product or in choosing an appropriate climbing wall, please contact us at wecare@funandfunction.com.



fun & function

www.funandfunction.com/safetyguidelines

Fun and Function LLC, PO Box 11, Merion Station, PA 19066

Safety Statements For Customer Assurance:

We earn your trust with our deep commitment and knowledge of the highest standards for children's product safety. In our design process, we optimize the features for safe use and functionality. Rigorous testing is conducted by independent and accredited third-party laboratories to meet or exceed government safety standards.

For manufacturing, we carefully select factories that share our dedication to product safety and quality control. Our Fun and Function team visits the factories and personally meets with the owners and workers -- and sometimes their families -- to ensure an understanding and respect for the children we serve.

We buy from the leading online retailer of high-quality products for children with sensory needs. We don't take shortcuts on quality or safety.

We do the hard work of finding good quality and useful products for you. We only work with suppliers with strong records of safety and commitment to quality. We stand behind our products. We don't take shortcuts with your children's safety.

We thoughtfully design our products to minimize the use of chemicals, like phthalates and BPA. Our team is dedicated to sourcing only high-quality products from trusted sources.

Our trained teams of occupational therapists and specialists help you find the right fit for your child's sensory needs.



fun & function

www.funandfunction.com/safetyguidelines

Fun and Function LLC, PO Box 11, Merion Station, PA 19066

Safety Statements:

Specific Categories for Website information and for customer assurance:

Bite Bands

Bite Bands are designed for mouthing. Tested by an independent and accredited third-party lab, this bite band meets or exceeds government safety standards for flammability and lead/heavy metals.

Chewies

Chewies are designed for mouthing. This chewy is made of FDA-approved silicone and is free of BPA. Tested by an independent and accredited third-party lab, this chewy meets or exceeds government safety standards for lead and phthalates. All chewy cords feature breakaway safety fasteners. While no chewy is indestructible, our chewies are tested to robust in-house standards of varying levels of pressure, for at least 2000 repetitive cycles, to meet the oral needs of the most active chewers. We provide use and care instructions, along with any applicable warnings, for all of our products in order to help our customers use the products properly and safely. Our professional team is always available to help with any questions you may have about installations, proper use and care, replacement parts, and returns before, during, and after your purchase.

Weighted Clothing

Weighted clothing is designed to support children with various sensory needs. Tested by an independent and accredited third-party lab, this clothing item meets or exceeds government safety standards for lead, flammability and physical and mechanical hazards. We also test the seams of the fabric and adjustable weights to help ensure the long-term durability for our most active wearers. We provide use and care instructions, along with any applicable warnings, for all of our products in order to help our customers use the products properly and safely. Our professional team is always available to help with any questions you may have about installations, proper use and care, replacement parts, and returns before, during, and after your purchase.

Weighted Blankets, Weighted Lap Pads and Weighted Companions

Weighted products are designed to support children with various sensory needs. Tested at an independent and accredited third-party lab, this product meets or exceeds government safety standards for lead, physical and mechanical hazards (and stuffing cleanliness). Seam strength is also tested to help ensure the long-term durability for our most active users. We provide use and care instructions, along with any applicable warnings, for all of our products in order to help our customers use the products properly and safely. Our professional team is always available to help with any questions you may have about installations, proper use and care, replacement parts, and returns before, during, and after your purchase.



fun & function

www.funandfunction.com/safetyguidelines

Fun and Function LLC, PO Box 11, Merion Station, PA 19066

Swings and Suspension Products

Swings, suspension products and hardware has been tested by an independent and accredited third-party lab and meets or exceeds government safety standards for physical and mechanical hazards, lead, load bearing and weight capacity. We provide use and care instructions, along with any applicable warnings, for all of our products in order to help our customers use the products properly and safely. Our professional team is always available to help with any questions you may have about installations, proper use and care, replacement parts, and returns before, during, and after your purchase.

Stretch fabric items

Stretch fabric products are designed to support children with various sensory needs. The fabric has been tested by an independent and accredited third-party lab and meets or exceeds government safety standards for flammability, physical and mechanical hazards. Where applicable, seam strength is tested for long-term durability for our most active users. We provide use and care instructions, along with any applicable warnings, for all of our products in order to help our customers use the products properly and safely. Our professional team is always available to help with any questions you may have about installations, proper use and care, replacement parts, and returns before, during, and after your purchase.

SensaSoft

This SensaSoft product, including the foam core and cover, has been tested by an independent and accredited third-party lab and meets or exceeds government safety standards for flammability, lead, phthalates, physical and mechanical hazards. We provide use and care instructions, along with any applicable warnings, for all of our products in order to help our customers use the products properly and safely. Our professional team is always available to help with any questions you may have about installations, proper use and care, replacement parts, and returns before, during, and after your purchase.

Gel Filled Products

This gel-filled product, including the outer and inner materials, has been tested by an independent and accredited third-party lab and meets or exceeds government safety standards for lead, phthalates, physical and mechanical properties, cleanliness, and flammability. The product is designed and manufactured for long-term durability and the gel liquid is non-toxic just in case the product wears down. We provide use and care instructions, along with any applicable warnings, for all of our products in order to help our customers use the products properly and safely. Our professional team is always available to help with any questions you may have about installations, proper use and care, replacement parts, and returns before, during, and after your purchase.

In a Pickle Game

This game has been tested by an independent and accredited third-party lab and meets or exceeds government safety standards for physical and mechanical hazards and lead.



fun & function

www.funandfunction.com/safetyguidelines

Fun and Function LLC, PO Box 11, Merion Station, PA 19066

Plush Products

Plush products have been tested by an independent and accredited third-party lab and meets or exceeds government safety standards for lead, seam strength, physical and mechanical hazards, flammability, small parts and stuffing cleanliness.

Scooters

This scooter has been tested by an independent and accredited third-party lab and meets or exceeds government safety standards for physical and mechanical hazards, lead, and phthalates. We provide use and care instructions, along with any applicable warnings, for all of our products in order to help our customers use the products properly and safely. Our professional team is always available to help with any questions you may have about installations, proper use and care, replacement parts, and returns before, during, and after your purchase.

Sensory Rooms/Sensory Room Products

All sensory room components have been tested by an independent and accredited third-party lab and meet or exceed government safety standards for flammability, physical and mechanical hazards, and lead. Electronic components are UL certified for safety. We provide installation, use, and care instructions, along with any applicable warnings, for all of our products in order to help our customers use the products properly and safely. Our professional team is always available to help with any questions you may have about installations, proper use and care, replacement parts, and returns before, during, and after your purchase.

Bubble Tubes, Tables and Walls

This product has been tested by an independent and accredited third-party lab, and meets or exceeds government safety standards for flammability, lead, and sharp points. Electronic components are UL certified for safety. We provide installation, use, and care instructions, along with any applicable warnings, for all of our products in order to help our customers use the products properly and safely. Our professional team is always available to help with any questions you may have about installations, proper use and care, replacement parts, and returns before, during, and after your purchase.

LED items

This product has been tested by an independent and accredited third-party lab and meets or exceeds government safety standards for flammability, physical and mechanical hazards, and lead. Electronic components are UL certified for safety.

Inflatables

This inflatable product has been tested by an independent and accredited third-party lab and meets or exceeds government safety standards for flammability, phthalates, lead, and physical and mechanical hazards. Our professional team is always available to help with any questions you may have about installations, proper use and care, replacement parts, and returns before, during, and after your purchase.



fun & function

www.funandfunction.com/safetyguidelines

Fun and Function LLC, PO Box 11, Merion Station, PA 19066

Flooring

This flooring has been tested by an independent and accredited third-party lab, and meets or exceeds government safety standards for flammability, lead and phthalates. Our professional team is always available to help with any questions you may have about installations, proper use and care, replacement parts, and returns before, during, and after your purchase.

Fidgets

This fidget product has been tested by an independent and accredited third-party lab and meets or exceeds government safety standards for lead, phthalates, physical and mechanical hazards, and small parts. Where appropriate, fidgets have been tested for seam strength. Our professional team is always available to help with any questions you may have about installations, proper use and care, replacement parts, and returns before, during, and after your purchase.

Sand Items

This weighted sand product has been tested by an independent and accredited third-party lab and meets or exceeds government safety standards for lead and physical and mechanical hazards. Sand products with small components are tested for small parts. The sand container has been tested for lead and phthalates requirements. Our professional team is always available to help with any questions you may have about installations, proper use and care, replacement parts, and returns before, during, and after your purchase.

Noise Reduction Headphones

These headphones have been tested by an independent and accredited third-party lab and meet or exceed government safety standards for lead content and for physical and mechanical hazards. Electronic components are UL certified for safety. Our professional team is always available to help with any questions you may have about installations, proper use and care, replacement parts, and returns before, during, and after your purchase.

Putty

This putty product has been tested by an independent and accredited third-party lab and meets or exceeds government safety standards for lead and physical and mechanical hazards. Putty products with small components are tested for small parts. The putty container has been tested for lead requirements. Our professional team is always available to help with any questions you may have about installations, proper use and care, replacement parts, and returns before, during, and after your purchase.

Cushions and Balls

Cushions and balls have been tested by an independent and accredited third-party lab, and meets or exceeds government safety standards for lead, physical and mechanical hazards, flammability and phthalates. Our professional team is always available to help with any questions you may have about installations, proper use and care, replacement parts, and returns before, during, and after your purchase.



fun & function

www.funandfunction.com/safetyguidelines

Fun and Function LLC, PO Box 11, Merion Station, PA 19066

Vibrating Products

This product, including a battery-operated vibration unit, has been tested by an independent and accredited third-party lab, and meets or exceeds government safety standards for lead, and physical and mechanical hazards. Electronic components are UL certified for safety. Our professional team is always available to help with any questions you may have about installations, proper use and care, replacement parts, and returns before, during, and after your purchase.

Foam Rollers

This product has been tested by an independent and accredited third-party lab and meets or exceeds government safety standards for lead and physical and mechanical hazards.

Scarves

This product has been tested by an independent and accredited third-party lab and meets or exceeds government safety standards for flammability.

Bean Bag Products

This bean bag product has been tested by an independent and accredited third-party lab and meets or exceeds government safety standards for flammability, lead, phthalates, physical and mechanical hazards and seam strength. We provide use and care instructions, along with any applicable warnings, for all of our products in order to help our customers use the products properly and safely. Our professional team is always available to help with any questions you may have about installations, proper use and care, replacement parts, and returns before, during, and after your purchase.

General statement

Our products have been tested by an independent and accredited third-party lab. Our products meet or exceed government safety standards for flammability, banned phthalates, lead, and physical and mechanical hazards. Our products comply with applicable voluntary, mandatory and best practice guidelines wherever such standards exist.



fun & function

www.funandfunction.com/safetyguidelines

Fun and Function LLC, PO Box 11, Merion Station, PA 19066